



Autodialer RLMC2009

The Advanced Autodial Unit can be connected to **Labcold** Monitoring systems and other laboratory equipment with suitable connections. The microprocessor based telephone dial out unit, once triggered by an external device, will dial up to four pre programmed telephone numbers and play a recorded warning message to alert you to any problems in your refrigeration equipment. It is made up of two components which can both be wall mounted. The Autodialer is installed by **Labcold** Service Technicians (or suitably qualified person).



Control unit

- Clear alphanumeric LCD display for status information and programming.
- Digital press-button keypad.
- Hinge down flap to prevent accidental keypad inputs.

Power Supply

- 240 volt AC input, 12-volt DC output.
- Integral self-recharging battery pack.
- Robust metal casing.
- 2m connecting cable to control unit with plug-in connections.

Not suitable for digital phone systems.

Part Number	RLMC2009
Configuration	Wall Mounted
Dimensions Control Module (H x W x D)	105 x 350 x 35mm
Dimensions Power Pack (H x W x D)	200 x 123 x 92mm
Connectivity	Direct, analogue line or extension telephone line
Power Connections	240 volt AC, 1 amp (Fused spur recommended)