


NEW Repeat Refrigerator Alarm from



Intelligently Designed for Safety and Security

Designed to work with the remote connections fitted as standard on the back of the Labcold IntelliCold® range of Pharmacy and Vaccine refrigerators, this battery operated repeat alarm is designed to allow you to hear and respond to a refrigerator alarm even if you are not in the same room. It is also ideal for noisy areas to boost the audible alarm fitted to the fridge.



Remote alarm contacts are located at the back of the fridge. The Repeat Refrigerator alarm kit is connected to the 2 pin connector on the left.

Supplied with a 5 meter long cable designed to work with the remote connections on the back of any refrigerator from the Labcold IntelliCold® range, the alarm is powered by a 9v battery, so no mains wiring is required.



www.labcold.com

