

SelfCookingCenter®



Connection set for three core temperature probes



Perfect results in iLevelControl.

Flexibility. But sure.

With iLevelControl, and especially with Finishing® with containers too, you can prepare different foods at the same time in your SelfCookingCenter®. In this way you make the most of your appliance and reduce your stress, even when there's a big rush. There is the "3 externally attached core temperature probes" option so you can now monitor the core temperature and the quality of your products even better. Thanks to accurate core temperature monitoring to the degree, you don't need to calibrate your goods anymore. Up to three different products are monitored at the same time. Regardless of whether you load the unit once or successively, you will always cook your food to the point. The option is available for all SelfCookingCenter® units with right hand doors, with the exception of the SelfCookingCenter® XS.

- 1 Rump steaks
core temperature **48°C**
- 2 Salmon fillets
core temperature **54°C**
- 3 Cutlets
core temperature **58°C**



Scope of delivery



Set with 2-point core temperature probes
For use in iLevelControl

Standard 6-point core temperature probe
For use in all cooking modes




Interactive display
Simply place the desired core temperature probe, select in the display and assign relevant tray - done. The probes are marked and also have different colours so that anyone can easily with them.



Connection socket outside the cooking chamber

Either connect the standard core temperature probe or the set of 3 core temperature probes. The socket can also be used for the RATIONAL sous-vide core temperature probe, available separately.

If no probe is inserted, the socket can be protected from spraying and dirt by the non-losable cover.



Holding box
When you are not using the core temperature probes, keep them safe in a holding box on the side of the unit which you can remove for cleaning and then re-attach.

Internationally tested



kiwa

DEKRA



RATIONAL International AG

Heinrich-Wild-Straße 202
9435 Heerbrugg
Switzerland

Tel. +41 (0)71 727 9090

Fax +41 (0)71 727 9080

info@rational-international.com

rational-international.com